Supporting Information

Cantilever Sensors for Rapid Optical Antimicrobial Sensitivity Testing

Isabel Bennett^{1,2}, Alice L. B. Pyne^{1,3*}, Rachel A. McKendry^{1,2*}

¹London Centre for Nanotechnology, University College London, 17-19 Gordon Street, London, WC1H 0AH UK ²Division of Medicine, University College London, Cruciform Building, Gower Street, London, WC1E 6BT UK ³Department of Materials Science and Engineering, Sir Robert Hadfield Building, University of Sheffield, Sheffield, S1 3JD, UK

*Corresponding authors: <u>r.a.mckendry@ucl.ac.uk</u>; a.l.pyne@sheffield.ac.uk

Pages: 9 Figures: 6 Tables: 1

SI Methods: Replication of Nanomechanical method.

Bacteria were cultured overnight at 37 °C, 250 r.p.m. ($\sim 10^6$ CFU/mL). The bacterial suspension was then washed in PBS, where 1 mL of bacterial suspension was centrifuged at 5000 r.p.m. for 1 minute. The supernatant was removed and the bacterial pellet re-suspended in 1 mL of PBS. This was repeated three times. The final bacterial solution was concentrated by a factor of 8x or 20x.

25% glutaraldehyde solution was diluted down to 0.5% in DI water. A small droplet was placed over cantilevers B and D (B: k = 0.12, $f_{res} = 23$ kHz; or D: k = 0.06, $f_{res} = 4$ kHz) of a DNP-S1 or NP-O10 chip (Bruker, USA) and incubated for 10 minutes. This was washed carefully with DI water and allowed to dry. A droplet of the concentrated bacterial solution was incubated on the same side for 30 minutes. Loose bacteria were washed off by dipping the cantilever gently into a petri dish of PBS. Bacterial immobilization cover was checked using a bright field microscope.

Experiments were carried out on a JPK NanowizardTM 3 ULTRA Speed AFM system (Bruker, USA) using DNP-S1 or NP-O10 cantilevers. The AFM was operated in contact mode for cantilever calibration. During experiments only the real-time scan function was used to monitor vertical deflection. Experiments were conducted at room temperature. Prior to the start of the experiments, the AFM laser was left on for ~12 hours to ensure the laser had warmed up fully and to reduce laser power fluctuations which would affect the drift of the signal.

Cantilevers were calibrated prior to bacterial immobilization. Sensitivity (nm/V) and spring constant (*k*) were measured by force curve analysis and thermal tune method, respectively. These were noted for conversion of raw mV data into nm.

The antibiotic, ampicillin, was used to kill the bacteria, and was added to a final concentration of 125 μ g/mL (far above the minimum bactericidal concentration (MBC)).

Supplementary Information Figures



Figure S1: Representative optical images of range in bacterial coverage. a, Low bacterial coverage. b, "Optimal" (500-600 cells) bacterial coverage. c, high 'clumpy' bacterial coverage. Scale bars = $50 \mu m$.



Figure S2: Investigation of bacterial immobilization conditions. Representative images of JM109 *E. coli*-immobilized cantilevers (DNP-S1 or tipless NP-O10), comparing different conditions and the number of bacteria immobilized. Other conditions are comparable across the cantilevers shown apart from the one specified. a, Cantilevers incubated with 8x (left) or 20x (right) concentrated bacterial solution for immobilization step. b, DH5 α (left) and JM109 (right) used as the bacteria for immobilization. c, Comparison of immobilization levels on the greater surface area of cantilever B versus the narrower cantilever D. d, Cantilevers treated with 0.5% glutaraldehyde (left) or 0.8% glutaraldehyde (right) prior to immobilizing bacteria. Bacterial immobilization number estimate shown bottom right of each figure. Scale bars = 50 µm. e, Box plot showing spread of data points across 28 immobilization experiments. IQR across these is below 500 cells, which is below "optimal" immobilization level for nanomechanical experiments. Median = 342 cells; Q₁-Q₃ = 238-531 cells.



Figure S3. Growth of bacteria over time is unaffected by addition of filtered media control. Number of bacterial crossings in filtered media (1), addition of more filtered media (2), inoculation with BL21-WT *E. coli* cells, and addition of more filtered media (4).



Figure S4. Baseline normalization and magnitude variability between experiments. **a**, Bacterial crossings over five experiments with BL21-WT *E. coli* (sensitive) and 125 μ g/mL ampicillin. Yellow, dark blue and orange experiments were inoculated with a higher CFU of bacteria than purple and light blue experiments. **b**, Normalized data for same experiments as **a**. Normalization of data to baseline pre-antibiotic treatment (blue highlighted section) to remove variability in bacterial inoculation level so experimental data comparable across experiments. **c**, r_{sensitivity} for five experiments shown in **a** and **b**. All show r_{sensitivity} < 1 (sensitive), demonstrating that initial inoculation level does not affect r_{sensitivity} value.



Figure S5. Kanamycin resistant and sensitive strain. Bacterial crossings over course of experiment for BL21-WT and BL21-kanaR + 125 μ g/mL kanamycin. r_{sensitivity} = 0.92 for BL21-WT (sensitive strain, green) and r_{sensitivity} = 2.0 for BL21-kanaR (resistant strain, red).



Figure S6. Data analysis steps applied to raw data. a, Raw data trace for vertical deflection taken over 800 seconds. **b**, Savitzky-Golay finite impulse response (FIR) smoothing filter of polynomial order 2 applied to the raw data, with a filtering frequency of 101 Hz. **c**, Peaks identified as having peak prominence value of 0.5 nm. This threshold was applied empirically across all files when carrying out the analysis to remove any bias of identifying peaks in the signal. **d**, Plot of bacterial crossings (number of peaks) across experiment. Each point is bacterial crossings collated for 267 seconds (800 seconds/3). r_{sensitivity} is calculated by taking ratio of S_{baseline} and 45 minutes post-antibiotic treatment, S_{antibiotic}. r_{sensitivity} provides a binary readout of sensitivity, r_{sensitivity} ≤ 1 indicates cell death or inhibition of bacterial growth, and sensitivity to the antibiotic in solution; r_{sensitivity} > 1 indicates bacterial growth, and therefore resistance to the antibiotic used. r_{sensitivity} for this example was 0.66.

| Antibiotic | Resistant/Sensitive |
|-------------------------|---------------------|
| Amikacin | R |
| Ampicillin | R |
| Augmentin | R |
| Cefotaxime | R |
| Ceftazidime | R |
| Chloramphenicol | R |
| Ciprofloxacin | R |
| Ertapenem | R |
| Fosfomycin | S |
| Gentamicin | R |
| Imipenem | S |
| Meropenem | R |
| Pipericillin/Tazobactam | R |
| Rifampicin | R |
| Temocillin | R |
| Tigecycline | S |
| Trimethroprim | R |
| Kanaymycin | R |

Table S1. Resistance spectrum of patient isolate from Great Ormond Street Hospital. Full resistance spectrum obtained using gold standard method for measuring break points for clinical isolate. Orange highlights two antibiotics chosen for this study, ampicillin and trimethoprim.